



**TITLE: EXAMINING IPTG-INDUCED  
EFFECTOR EXPRESSION LEVELS IN  
EPEC COMPLEMENTED STRAINS BY  
WESTERN BLOTTING**

Date: 12.10.2021

By: Ipsita Nandi & Efrat Zlotkin-Rivkin

1. Streak a plate for EPEC strain containing the desired effector protein encoding plasmid (according to the EPEC strains inventory file) with appropriate antibiotics and incubate O/N @ 37°C. **IMPORTANTLY:** These EPEC strains are typically mutants deleted in a specific effector, complemented with a plasmid encoding an effector that is typically tagged (HA, FLAG, SBP etc) – e.g. EPEC-*espH* + EspH-SBP. Therefore, the control experiment for these experiments will be the EPEC bearing the deleted gene; e.g. EPEC-*espH*.
2. **LB starter culture:** Pick a single colony from the plate in step-1 and inoculate into a snap cap tube with 3mL LB+ appropriate antibiotics. **Incubate O/N in 37°C without shaking.** Include a tube with LB and no bacteria as control for unwanted infection- **if this control is turbid after O/N culturing, then your LB is considered contaminated- discard the experiment!**
3. **Activation tube:** Aliquot 5 ml of an appropriate tissue culture medium (usually HDMEM) in 50 ml centrifuge tube, the cap should be loosely closed and place it in CO<sub>2</sub> incubator for overnight. **Media should not contain any antibiotics and serum.**
4. Following day, inoculate 100 µl of the **LB starter culture** into the **activation tube** containing 5 mL of the appropriate tissue culture medium (1:50 ratio; LB starter culture: Tissue culture media, typically DMEM or MEM).
5. Incubate the culture from **Step 4** in 37°C + 5% CO<sub>2</sub> (without shaking) for 3 hours. Pay attention to keep the test tube cap loose to allow the entry of CO<sub>2</sub>; After 2.5 hrs. add the appropriate concentration of IPTG and place it back into the incubator for another half an hour. Use **freshly prepared** IPTG stock diluted to **0, 0.05, 0.1, 0.2, 0.25, 0.5** mM into the DMEM/activated bacteria mixture. There is no need to add IPTG to the deletion mutant.

6. After the incubation is over, mix the tube and immediately take 1 ml of the culture into 1.5 ml micro-centrifuge tube, spin at 6000 rpm for 2 mins.
7. Discard the supernatant and resuspend the pellet in 200  $\mu$ l 1X SDS-PAGE sample buffer.
8. Heat 95 °C for 10 min, spin-down maximum speed for 1 min.
9. Take sup and load 20  $\mu$ l per sample well onto a gel (since the effectors are small – up to 20 kDa – 13% gels are typically used) and analyze by Western blotting.

Gels should be loaded in the following order:

1. Mw markers
2. EPEC-*espH*
3. EPEC-*espH*+EspH-SBP (0 mM IPTG)
4. EPEC-*espH*+EspH-SBP (0.05 mM IPTG)
5. EPEC-*espH*+EspH-SBP (0.1 mM IPTG)
6. EPEC-*espH*+EspH-SBP (0.2 mM IPTG) or (0.25 mM IPTG)
7. EPEC-*espH*+EspH-SBP (0.4 mM IPTG) or (0.5 mM IPTG)