

TITLE: INVASION ASSAY:

Name: Rachana Ramachandran

Experiment:

- 1. Culture HeLa cells in DMEM supplemented with 10% FCS and 1% PSA and glutamine at 37 °C in a CO₂ incubator.
- 2. Count and seed appropriate number of cells in 12-well plate (40,000 cells/well for 2 day and 80,000 cells for 1 day).
- 3. [In KD experiments (with ABR siRNA), 24-well plate was used (12,000 cells/well was added and after transfection with ABR siRNA, cells were kept for 72 hours at 37 °C in a CO₂ incubator.
- 4. Start overnight culture of bacteria in 3 ml of LB broth with appropriate antibiotics and incubate at 37 °C with no shaking.
- 5. Add bacteria from the overnight culture (at a ratio of 1:50) to warm plain DMEM and incubate for 3 hours at 37 °C in a CO₂ incubator.
- 6. Wash cells 2x with 1xPBS
- 7. Add 1ml of activated bacteria to each well and incubate for 90 min at 37 °C in a CO₂ incubator.
- 8. Remove the medium and add instead 1 ml DMEM supplemented with 100 μg/ml gentamicin (Cat#G1264, Sigma Aldrich) or without gentamicin (as control).
- 9. Incubate for 90 min at 37 °C in a CO₂ incubator.
- 10. Wash cells 3x with 1xPBS.
- 11. Lyse the cells with 1 ml lysis buffer containing 10 mM Tris (pH-7.4) and 1% TritonX100.
- 12. Mix the lysates well by few rounds of up and down pipetting.
- 13. Serially dilute the lysate to few different dilutions in 1xPBS (1:500; 1:1000; 1:2000; 1:4000; 1:8000).
- 14. From each dilution, $100\mu l$ were plated in LB-agar plates and incubated overnight at $37^{\circ}C$.
- 15. Count the number of colonies in each plate with colony counter.

Quantification:

- 1. The concentration of protein in each sample is calculated by BCA assay kit method.
- 2. Number of colonies in each plate for cells treated with gentamycin (invaded bacteria) were counted and multiplied to its respective dilution to determine the cfu/ml (colony forming unit/ml). The number of colonies is divided with its protein concentration.
- 3. The number of colonies in all the plates, for cells not treated with gentamycin (total bacteria) serves as control (were equal in all cases, 8×10^6 cfu/ml).